

Superior Irrigation and Evacuation for Endoscopy



Single use tubing with unique geometry head.



Control Bowel Poor-Prep. Clear Lower GI Bleeding.

Fast, convenient, reliable, clear view

Challenge: To Proceed Every Time

Colonoscopy preparation is uncomfortable and time consuming. More than ten percent of patients arrive for their procedure under-prepared ... at increased risk of false negative or aborted procedure.

Solution: Use ClearPath Lower GI

ClearPath Lower GI single use tubing with unique geometry head is cost-effective every colonoscopy, to reliably control bowel poor-prep.

- Complete colonoscopy easily and cost effectively, even if patient is poorly prepped.
- Reduce poor-prep aborted exams to near-zero ... help protect doctor-patient relationship.
- Minimize patient discomfort and healthcare expense associated with aborted procedure.
- Decrease false negatives by overcoming suboptimal visibility.
- Keep colonoscope working channel free for instrumentation.
- Easily pair with all major colonoscopes.
- No need to change procedural technique.



Potent irrigation dislodges debris.



Potent evacuation clears visual field.



Superior Irrigation and Evacuation for Endoscopy

Control Bowel Poor-Prep. Clear Lower GI Bleeding.


Problem: Identify Source and Control Bleeding

Lower GI bleeding calls for high-volume evacuation to remove blood and large clots, along with working-channel access to the identified source.


Solution: Use ClearPath Lower GI

ClearPath Lower GI provides large cross section evacuation while keeping colonoscope working channel available every moment of procedure ... exploit clear view with intervention to control bleeding.

- Identify source and control lower GI bleeding.
- Keep colonoscope working channel free for instrumentation.
- Easily pair with all major colonoscopes.
- No need to change procedural technique.



Potent irrigation displaces blood and large clots.



Potent evacuation clears visual field.

Product Specification	
Irrigation pressure	1.5 - 2.5 bar
Evacuation rate	1800 ml/min
Evacuation cross section	27 mm ²
Use with	ClearPath Controller

EasyGlide Ltd.






30 Ha'Ella St., Kfar Truman 73150, Israel Fax/Phone: 972.77.6201003
info@easy-glide.com • www.easy-glide.com

ClearPath is a trademark of Easy-Glide Ltd.

©2011 Easy-Glide, Ltd. All rights reserved