

EasyGlide Product List

Product	Type	Catalog number	Tip diameter	Quantity per box	Notes
ClearPath Upper GI	Small Upper GI	CPU 4800	8.8-9.8 mm	20	
	Medium Upper GI	CPU 4700	9.8-10.8 mm	20	
ClearPath Lower GI	Small (Pediatric) - Long	CPP 1600	11.0-12.0 mm	20	1680 mm length
	Small (Pediatric) - Short	CPP 1610	11.0-12.0 mm	20	1330 mm length
	Medium (Adult) - Long	CPA 1500	12.7-13.4 mm	20	1680 mm length
	Medium (Adult) - Short	CPA 1510	12.7-13.4 mm	20	1330 mm length
	Large (Adult) - Long	CPA 1700	13.5-14.2 mm	20	1680 mm length
	Large (Adult) - Short	CPA 1710	13.5-14.2 mm	20	1330 mm length
ClearPath Tubing (for Irrigation)	---	CPT 2200	All scopes	20	
ClearPath Irrigator	Green (Olympus/Fujinon) Adaptor	CPI 5100	Olympus, Fuji	50	(see matrix for supported models)
	Purple (Pentax) Adaptor	CPI 5200	Pentax	50	
ClearPath Controller	---	CC 3100		1	
ClearPath Pole Cart	---	ET 6000		1	

The above products are for sale outside of the USA. For the list of products cleared for marketing in the USA, please contact a sales representative.



30 Ha'Ella St., Kfar Truman 73150, Israel Fax/Phone: 972.77.6201003
info@easy-glide.com • www.easy-glide.com

ClearPath is a trademark of Easy-Glide Ltd.

©2011 Easy-Glide, Ltd. All rights reserved

Lower GI Scope/Device Table

ClearPath device	Scope	Scope size
Small (Pediatric) Lower GI – Short CPP 1610	Fuji EC-530MP	11.0mm/1330mm
	Olympus PCF-Q180A / H180A / Q160A (I)	11.3mm/1330mm
	Pentax EC-3490 Mi	11.6mm/1300mm
Small (Pediatric) Lower GI – Long CPP 1600	Fuji EC-530LP	11.0mm/1690mm
	Olympus PCF-Q180A / H180A / Q160A (L)	11.3mm/1680mm
	Pentax EC-3490 / 3485 / 3430-LK / Li / TLi Pentax EC-3490 -FK / Fi	11.6mm/1700mm 11.6mm/1500mm
Medium (Adult) Lower GI – Short CPA 1510	Olympus CF-Q165 / Q180A / Q160 / Q160A (I) Olympus CF-Q145 (I)	12.8mm/1330mm
Medium (Adult) Lower GI – Long CPA 1500	Olympus CF-Q165 / Q180A / Q160 / Q160A (L) Olympus CF-Q145 (L) / CF-Q260A (L)	12.8mm/1680mm
	Fuji EC-530WI Pentax EC-3830LK	12.8mm/1520mm 12.8mm/1700mm
Large (Adult) Lower GI – Short CPA1710	Olympus CF-H180A / Q / D160D (I)	13.8mm/1330mm
Large (Adult) Lower GI – Long CPA 1700	Olympus CF-H180A / Q / D160D (L)	13.8mm/1680mm

 EasyGlide Ltd.

 ClearPath™ Controller
  ClearPath™ Upper GI
  ClearPath™ Irrigator
  ClearPath™ Tubing

30 Ha'Ella St., Kfar Truman 73150, Israel Fax/Phone: 972.77.6201003
info@easy-glide.com • www.easy-glide.com

ClearPath is a trademark of Easy-Glide Ltd.

©2011 Easy-Glide, Ltd. All rights reserved

Upper GI Scope/Device Table

ClearPath device	Scope	Scope size
Small Upper GI CPU 4800	Olympus GIF-Q180	8.8mm/1030mm
	Pentax EG-2790K	9.0mm/1030mm
	Olympus GIF-Q165	9.2mm/1030mm
	Fuji EG-250WR5, EG-250CT5, EG-530WR	9.4mm/1030mm
	Olympus GIF-Q160 / H260 / Q260	9.5mm/1030mm
Medium Upper GI CPU 4700	Olympus GIF-H180 / H180J / Q140 / Q260J	9.8mm/1030mm
	Pentax EG-201 FP	9.8mm/1100mm
	Pentax EG-2990K / 2985K / 2930K	9.8mm/1050mm

Important note and disclaimer:

- Inconsistencies may occur due to changes or variations in the actual scope size
- The tip size of refurbished scopes may be different than the value specified in the table and the product documentation

Irrigator GI Scope/Device Table

ClearPath device	Scope
CPI 5100 Olympus and Fujinon (Green)	Olympus
	Fuji (except some models such as EG-201 FP, EC-201 WL)
CPI 5200 Pentax (Purple)	Pentax

 **EasyGlide** Ltd.






30 Ha'Ella St., Kfar Truman 73150, Israel Fax/Phone: 972.77.6201003
info@easy-glide.com • www.easy-glide.com

ClearPath is a trademark of Easy-Glide Ltd.

©2011 Easy-Glide, Ltd. All rights reserved